

Clinical cases by  
Dr. Fernando Rojas-Vizcaya


dental  
bone & tissue  
regeneration




botiss  
biomaterials

*strictly biologic*

# botiss BTR system: BONE



# botiss BTR system: TISSUE



# Dr. Fernando Rojas-Vizcaya, DDS, MS

Castellon, Spain

Dr. Fernando Rojas-Vizcaya is Adjunct Assistant Professor in the Department of Prosthodontics at the University of North Carolina in Chapel Hill, NC. He is the Founder and Director of both; the Mediterranean Prosthodontic Institute ([www.prosthodontics.es](http://www.prosthodontics.es)) and BoneModels ([www.bonemodels.es](http://www.bonemodels.es)), and maintains a private practice limited to dental implant surgery and prosthodontics in Castellon, Spain, and he lectures extensively worldwide. His major research interest includes esthetic management in complex dental implant cases in immediate placement and immediate loading protocols.



## Education

Post-Doctoral Implant Research Fellow	2005-2006	University of North Carolina, USA
Specialty Degree in Prosthodontics	2002-2005	University of North Carolina, USA
Master of Science in Prosthodontics	2002-2004	University of North Carolina, USA
Doctorate in Medicine and Bucco Facial Surgery	1998-2000	University Complutense of Madrid, Spain
Specialty in Osseointegrated Implants	1994-1995	University Complutense of Madrid, Spain
Specialty in Oral Surgery	1993-1994	Hospital Gregorio Marañón of Madrid, Spain
Specialty in Oral Medicine	1992-1993	University Complutense of Madrid, Spain
DDS Degree	1985-1989	University Javeriana of Bogota, Colombia





# Clinical cases

Implantology, Periodontology, Oral Surgery & CMF

Extraction socket, vertical augmentation of the bone defect with maxresorb<sup>®</sup> inject and Jason<sup>®</sup> membrane



Clinical situation before augmentation with vertical bone loss



Visible bone resorption following extraction of canine



Clear bone defect visible on buccal bone site



Augmentation with maxresorb<sup>®</sup> inject



Bone defect is filled with maxresorb<sup>®</sup> inject



maxresorb<sup>®</sup> inject covers the buccal wall bone defect



Buccal wall protected with Jason<sup>®</sup> membrane



Jason<sup>®</sup> membrane turned down over the socket




Final suturing of the socket

# Clinical cases

Implantology, Periodontology, Oral Surgery & CMF

Extraction socket with bone defect; GBR with Jason<sup>®</sup> membrane and cerabone<sup>®</sup>



Lateral CT scan of the mandible showing the buccal perforation in the middle of the root of the lateral and the bone defect



Clinical situation before tooth extraction



Clinical situation after tooth extraction



Bony defect of the buccal bone wall visible after soft tissue mobilization



Buccal bone wall protected with Jason<sup>®</sup> membrane



Socket is filled with cerabone<sup>®</sup> particles



Jason<sup>®</sup> membrane turned down over the socket



Wound closure and suturing



Wound closure and suturing



## Clinical cases Implantology, Periodontology, Oral Surgery & CMF

Extraction socket, fenestration defect, ridge augmentation;  
maxgraft®, mucoderm® and Jason® membrane,



Situation before extraction of the teeth



Situation before extraction of the teeth, occlusal view



Atraumatic tooth extraction



Fenestration defect visible after extraction of the teeth



Situation with extracted teeth before augmentation



Rehydration of mucoderm® matrix and maxgraft® allogenic bone



Placement of trimmed mucoderm® for soft tissue augmentation




Jason® membrane cut to shape for vestibular placement




Covering of the vestibular wall with Jason® membrane

## Clinical cases Implantology, Periodontology, Oral Surgery & CMF



Jason® membrane protecting the vestibular wall



Extraction sockets filled with maxgraft® granula




Jason® membrane turned down over the augmented area



Jason® membrane covering the augmented area



Fixation of mucoderm® and Jason® membrane by sutures



Fixation of mucoderm® and Jason® membrane, occlusal view



Final suturing of the sockets



Final suturing of the sockets; occlusal view



Provisional removable partial denture without compression



Clinical cases Implantology, Periodontology, Oral Surgery & CMF

Extraction socket; collacone®



Situation before extraction



Situation following atraumatic tooth extraction, alveolar walls intact



Filling of socket with collacone®



Suturing and open healing



Clinical situation after healing



No bone resorption visibel after opening of the flap



Placement of implant after 3 months



Placement of healing abutment and suturing



Clinical situation after implantation and healing after 2 months

Clinical cases Implantology, Periodontology, Oral Surgery & CMF

Socket preservation and augmentation with cerabone® and Jason® membrane



Maxilla showing a fixed partial restoration



Removal of bridges



Situation after removal of bridges



Bridges and tooth roots extracted from the alveolae



Alveoles after extraction of tooth roots




Covering of the lingual wall with Jason® membrane



Filling of sockets and augmentation with cerabone®



Covering of the augmentation area with Jason® membrane



Jason® membrane turned down over the augmentation area



Suturing of the membrane



Wound closure and suturing



Clinical situation after 2 weeks healing




# Clinical cases

Implantology, Periodontology, Oral Surgery & CMF

Extraction socket and fenestration defect;  
Jason® membrane, Jason® fleece and cerabone®



Situation before extraction



DVD scan of the maxilla showing fenestration defect



Atraumatic tooth extraction




Atraumatic tooth extraction



Clinical situation after tooth extraction



Clear fenestration defect visible at the lingual side



Jason® membrane covering the bone defect



Block of the apex of the socket with Jason® fleece



Filling of the sockets with cerabone®



Covering of sockets with Jason® fleece



Final suturing of the sockets



Provisional removable prosthesis without compression in the area of the sockets

# Clinical cases

Implantology, Periodontology, Oral Surgery & CMF

Socket preservation with Jason® membrane, cerabone® and Jason® fleece



Clinical situation before tooth extraction and implantation



Clinical situation after drilling of implant bed



Situation after implantation and tooth extraction, visible bone defect



cerabone® and Jason® fleece placed in the extraction socket



Covering of the buccal bone defect with Jason® membrane



Final wound closure and suturing

Socket preservation with Jason® membrane and cerabone®



Clinical view before surgery



Buccal bone defect visible after soft tissue mobilization



Placement of implants



Covering of the defect area augmented with cerabone® with Jason® membrane



Perfect adaptation of Jason® membrane to the bone



Immediate provisionalization using prefabricated acrylic crowns



botiss dental GmbH  
Knesebeckstrasse 59-61  
10719 Berlin / Germany  
Fon +49 30 20 60 73 98 30  
Fax +49 30 20 60 73 98 20  
[contact@botiss.com](mailto:contact@botiss.com)  
[www.botiss.com](http://www.botiss.com)  
facebook: botiss biomaterials



botiss  
biomaterials